

FREIBURG

Chair weigher up to 150 kg body weight

Equipment:

- aluminium back with PUR backrest
- swinging armrests enable easy side transfer
- upholstered seat
- swinging foot support
- swivelling castors, in the rear with parking/direction brake
- approved up to 150 kg body weight

Indicator:

- backlighted LCD indication
- auto-off function (time adjustable)
- units in kg or lb
- weighing range 150 kg
- division 100 g

Types:


- ref. no. 349.30.01 battery or plug operation
- ref. no. 349.30.02 accumulator operation incl. charger
- Art.-Nr. 349.30.03 calibrated - battery or plug operation
- Art.-Nr. 349.30.04 calibrated - accumulator operation incl. charger

30 grey/fire-red

Up to 150 kg body weight!



ref. no. 349.30.01



protected by patent

Technical data:

Type:	Freiburg
depth:	95 cm
width:	63 cm
seath width (A):	50 cm
seat height (B):	55 cm
seat depth (C):	45 cm

Type:	Freiburg
weight:	18 kg
max. body weight:	150 kg
division:	100 g
ref. no.:	349.30. ..

Technical changes and changes in colour reserved.

Chair weigher